

SEROLOGICAL WATER BATH

PRODUCT DESCRIPTION

WE ARE A RENOWNED ENTITY IN THIS DOMAIN, INVOLVED IN MANUFACTURING, EXPORTING AND SUPPLYING SEROLOGICAL WATER BATH. THIS BATH IS MANUFACTURED UNDER THE STRICT DIRECTION OF OUR EXPERIENCED PROFESSIONALS USING EXCEPTIONAL QUALITY COMPONENTS WITH THE AID OF CONTEMPORARY TECHNIQUES. THE OFFERED BATH IS COMMONLY USED IN REGULAR BACTERIOLOGICAL AND LABORATORY APPLICATIONS REQUIRING INCUBATION AND GENERAL TEST PROCEDURES. MOREOVER, OUR VALUED CLIENTS CAN AVAIL THIS SEROLOGICAL WATER BATH IN VARIOUS DIMENSION AND CAPACITY AT MARKET LEADING PRICES.

FEATURES:

- LONG LIFE
- TROUBLE FREE PERFORMANCE
- USER FRIENDLY DESIGN
- OPTIMUM FUNCTIONALITY

SPECIFICATIONS



MODELS	INNER DIMENSIONS (L X W X D) APPROX.	CAPACITY APPROX.	POWER	NO. OF TUBES
Cow -T	300 X 250 X 175 MM	13 LTRS.	1.5 KW	2
Cow -T	375 X 300 X 175 MM	18 LTRS.	2.0 KW	4
Cow -T	455 X 300 X 175 MM	24 LTRS.	2.0 KW	6
Cow -32T	605 X 300 X 175 MM	32 LTRS.	4.0 KW	8

Rectangular Water Bath

The electrically operated Rectangular Water Bath are made of thick stainless steel and feature constant water levels maintained inside them. For optimum performance, these also come complete with concentric rings having 75 mm dia, wire plugs and with suitable drainage facility. The bath machine also comprises thermostatic controller that allows better control of temperature levels with removable glass tray and stainless steel platform for drying of sections.

Features:

Versatile

Coated with an anti resistant material

Are able to give lifelong performance



Specification:

6 holes of 75 mm dia.

12 holes of 75 mm dia.